

Polarographic Analyses Half wave potentials of metal ions

Of interest to:
General analytical chemistry
B 1

Summary

In the following tables the half-wave potentials or peak potentials of 90 metal ions are listed. The half-wave potentials (listed in Volt) are measured at the dropping mercury electrode (DME) at 25°C unless indicated otherwise.

Electrodes

WE	MME	
AE	Pt or Glassy Carbon	6.1246.020 6.0343.000 6.1241.020 + 6.1247.000
RE	Ag/AgCl/KCl (3 mol/L)	6.0728.020 6.1245.010

In some cases other electrodes are specified.

Legend

- w Well marked wave or peak
- i Poor wave or evaluation impossible (e.g. in the region of increasing background current of the supporting electrolyte)
- cw Catalytic reduction of hydrogen (in some cases suitable for quantitative analysis)
- NR No signal of the respective metal ion (polarographic or voltammetric analysis not possible)
- ↓ Formation of precipitate. Attention to solubility product!
- No information available

Part 1: Various metals and transition metals

Electrolytes

Acetate	c(CH ₃ COOH) = 2 mol/L + c(CH ₃ COONH ₄) = 2 mol/L pH 4.6
NH ₃ / NH ₄ Cl	c(NH ₃) = 1 mol/L + c(NH ₄ Cl) = 1 mol/L pH 9.6
KCl	c(KCl) = 1 mol/L
HCl	c(HCl) = 1 mol/L
NaOH	c(NaOH) = 1 mol/L
KNO ₃	c(KNO ₃) = 0.1 mol/L
KSCN	c(KSCN) = 1 mol/L
TEA	c(Triethanolamine) = 0.3 mol/L + c(KOH) = 0.1 mol/L
Tartrate	c(Na ₂ -tartrate) = 0.5 mol/L pH = 9.0
Na ₂ EDTA	c(Na ₂ EDTA) = 0.1 mol/L

Citrate	c(Na ₃ -citrate) = 0.1 mol/L + c(NaOH) = 0.1 mol/L
Oxalate pH 2	c(Ammonium oxalate) = 0.1 mol/L, pH setting with H ₂ SO ₄
Oxalate pH 4	c(Ammonium oxalate) = 0.1 mol/L, pH setting with H ₂ SO ₄
Oxalate pH 6	c(Ammonium oxalate) = 0.1 mol/L, pH setting with H ₂ SO ₄
Oxalate pH 8	c(Ammonium oxalate) = 0.1 mol/L, pH setting with NH ₃

Metal Ion	Acetate	NH ₃ /NH ₄ Cl	KCl	HCl	NaOH	KNO ₃	KSCN	TEA	Tartrate	Na ₂ EDTA	Citrate	Oxalate pH 2	Oxalate pH 4	Oxalate pH 6	Oxalate pH 8
Al(III)	NR	NR	-1.59 i	NR	NR	NR	NR	NR	NR	NR	NR	-	-	-	-
As(III)	-0.95	-1.62 i	-0.07 i	-0.38 w -0.62	-0.25	-0.11 w	NR	-0.16	NR	-1.23 w	-0.22 w	NR	-1.56 w	-1.70 i	-1.58 i
Bi(III)	-0.21 w	NR	-0.06 w	-0.07 w	-0.68	-1.17 i ↓	NR	-0.87 w	-0.31 w cw	-0.54 w	-0.77 w	-0.13 w	-0.21 w	-0.22 w	-0.26 w ↓
Cd(II)	-0.63 w	-0.77 w	-0.62 w	-0.61 w	-0.82 i	-0.55 w	-0.62 w	-0.81	-0.64	NR cw	-0.68 i	-0.56 w	-0.61 w	-0.62 w	-0.62 w
Ce(IV)	NR	NR	NR	NR	NR	NR ↓	NR	NR	NR	NR	NR	-	-	-	-
Co(II)	NR	-1.26	-1.38	NR	NR ↓	-1.27 w -1.54	-1.11	NR	-1.53	NR	NR	NR	NR	NR	-1.52 i
Cr(III)	-1.20 i	-1.34 w ↓	-1.03	-1.00 w	NR ↓	-0.87 -1.01	-1.01 w	NR	-1.58 i	-1.20 w	NR	NR	NR	-1.73 i	NR
Cr(VI)	NR	-0.28 w -1.45 w -1.68	-0.38 -0.88 w	-1.03 i	-0.85 w	-0.32 -1.10	-0.41 -1.00 w	-1.27 w	-0.33 -0.86 -1.71	-1.20 w	-0.73 w	NR	NR	-0.02 w	-1.14 w -1.72 i
Cu(II)	-0.25 w	-0.20 w -0.46 w	-0.18	-0.20	-0.35 -0.43	NR	-0.64	-0.50	-0.08 w	-0.26 w	-0.41 w	-0.05 w	-0.15 w	-0.18 w	-0.18 w
Fe(II)	NR	-1.44 w	-1.57 i cw	i cw	-1.55 i ↓	-1.31 -1.56 i cw	-1.59 w	-1.02 w -1.68	-1.51 i	-0.09 w	-0.83 w -1.57 w	-0.02 w	-0.17	-0.18 w	-0.18 w -1.61 i
Fe(III)	-0.02 -0.24	NR ↓	-1.42 i	NR	NR ↓	-1.29 i	-1.56 w	-0.99 w -1.68	-0.19 -1.53 w	-0.11 w	-0.85 w -1.59 w	-0.18 w	-0.17 w	-0.19 w	-0.19 w
Ge(IV)	-1.34 i	-1.41	NR	NR	NR	-1.26 i cw	i cw	NR	NR	NR	NR	NR	NR	-1.38 w	-1.38 w
In(III)	-0.67 w	-1.08	-0.57	-0.57	-1.12 w	-0.90 i	-0.59 w	-1.51 w	NR	NR	-1.17	-0.61 w	-0.67 w -0.84 w	-0.68 w -0.87 w	-0.70 w -1.18
Mn(II)	NR	-1.59 w	-1.55 w	NR	-1.68 i ↓	-1.46 i	-1.56 i	-0.44 -1.66 i	-1.52 w	NR	NR	NR	NR	-1.61 i	-1.60 i

Metal Ion	Acetate	NH ₃ /NH ₄ Cl	KCl	HCl	NaOH	KNO ₃	KSCN	TEA	Tartrate	Na ₂ EDTA	Citrate	Oxalate pH 2	Oxalate pH 4	Oxalate pH 6	Oxalate pH 8
Mo(VI)	-0.63 w -1.15 -1.30	NR	NR	-0.08 i	NR	NR	NR	NR	NR	-0.52 -0.75 w	NR	-0.02 -0.11 -0.33 w -0.74 w	-0.17 w -0.29 -0.59	NR	NR
Ni(II)	-1.09 w	-1.07 w	-1.03	NR	NR ↓	-1.00 w	-0.66 w	-1.36	NR	NR	NR	NR	NR	NR	-1.32
Pb(II)	-0.47 w	-0.47 w	-0.41 w	-0.41 w	-0.73 ↓	-0.36 w	-0.39 w	-0.91	-0.59	-1.03 w	-0.69 w	-0.40 w	-0.45 w	-0.56 w	-0.57 w
Pd(II)	-0.53 -0.76	-0.75 w	i cw	NR	-1.33 i	NR	i cw	-1.01 i ↓	NR	NR	NR	-	-	-	-
Sb(III)	-0.39 -0.50 w	-0.81 w	-0.07 -0.15 w	-0.12 w	-0.42 -1.18	-0.19 w	-0.58 w -0.70	-1.31 w	NR	-0.64 w	-0.33 -0.95 w	-0.33 w	-0.46 w	-0.41 w -0.57 w	-0.76 w
Se(IV)	-0.70 w -1.19	-1.54 w	NR	-0.06 -0.47	NR	NR	NR	NR	-1.21	-0.67 w -1.20 w	NR	-0.03 -0.60 w	-0.68 w -1.11	-1.35	-1.35
Sn(II)	-0.13 w -0.63 w	-0.72	-0.43 w	-0.45 w	-0.83 -1.15 w	-0.35 w	-0.45 w -1.59 w	-0.75 w -1.15 w	-0.54 -0.83 -1.06	-0.13	-0.83 w -1.08 w	-0.12 w -0.51 w	-0.15 w -0.62 w	-0.26 w -0.62 w	-0.51 i -0.64 i
Sn(IV)	NR	NR	NR ↓	-0.45	NR ↓	NR	-1.57	NR	NR	-1.16 i	NR	-0.51 w	NR	NR	NR
Ti(III)	-	-	-	-	-	-	-	-	-	-	-	-0.31 w	-0.41 w	-0.52	-0.70
Ti(IV)	-	-	-	-	-	-	-	-	-	-	-	-0.29	-0.42	NR	NR ↓
Tl(I)	-0.44 w	-0.46 w	-0.47	-0.46 w	-0.45	-0.43 w	-0.50 w	-0.45 w	-0.45 w	-0.44 w	-0.47 w	-0.44	-0.44	-0.44	-0.44
U(VI)	-0.42 w	NR	-0.18 i	-0.17	-0.92 i	NR	NR	-0.87	NR	NR	NR	-0.18 w	-0.32	-0.38 w	-0.42 w
V(V)	i cw	-1.12 -1.29 i	-1.10 i cw	-1.11 i	-0.40 i	-1.03 i	-0.50 w	-0.35	-0.33 i	-1.23	-0.76 -1.11 w	NR	NR	-0.90 -1.30 w	-0.19 w -0.88 -1.29 w
W(VI)	-0.26 i	NR	NR	-1.01 i cw	NR	NR	NR	NR	NR	-1.26 i	NR	-	-	-	-
Zn(II)	-1.05	-1.32 w	-0.99 w	-0.99 i	-1.58	-0.97	-1.01 w	-1.58	-1.27 w	NR	-1.38 w	NR	NR	-1.38	-1.29
Zr(IV)	NR	NR	-1.49 w	NR	NR	-0.95 w	-1.58 w	NR	NR cw	NR	NR	-	-	-	-

Part 2: Alkaline and alkaline earth metals

Electrolytes

c(TEAI) = 0.1 mol/L	c(tetraethylammonium iodide) = 0.1 mol/L in water
c(TEAI) = 0.1 mol/L, 20% DMF	c(tetraethylammonium iodide) = 0.1 mol/L in water:dimethylformamide 80:20

Half-wave potentials

Metal ion	Supporting electrolyte	Half-wave potential
Ba(II)	c(TEAI) = 0.1 mol/L	-1.90
Be(II)	c(TEAI) = 0.1 mol/L, pH 2.3	-1.91
Ca(II)	c(TEAI) = 0.1 mol/L	-2.24 i
Cs(I)	c(TEAI) = 0.1 mol/L	-2.04
K(I)	c(TEAI) = 0.1 mol/L	-2.11
Li(I)	c(TEAI) = 0.1 mol/L	-2.31
Mg(II) ¹	c(TEAI) = 0.1 mol/L	²
Na(I)	c(TEAI) = 0.1 mol/L	-2.07
Rb(I)	c(TEAI) = 0.1 mol/L	-2.06
Sr(II)	c(TEAI) = 0.1 mol/L, 20% DMF	-1.85

Remarks

Alkaline and alkaline earth metals usually cannot be determined in real samples because the signals are not well separated from each other.

¹ Magnesium is determined indirectly as an oxime or as a complex with solochrome violet RS.

² Maximum, non exploitable.

Part 3: Rare-earth elements

Electrolytes

c(TEA-CIO ₄) = 0.1 mol/L in DMF	c(tetraethylammonium perchlorate) = 0.1 mol/L in dimethylformamide
---	--

Half-wave potentials

Metal Ion	Supporting electrolyte	Half-wave potential ³
Ce	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.44
Dy	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.08
Er	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.16
Eu	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-0.80 -2.30
Gd	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.12
Ho	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.20
La	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.16
Lu	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.08
Nd	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.22
Pm	c(TEA-CIO ₄) = 0.1 mol/L in DMF	– (radioactive)
Pr	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.24
Sm	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-1.98
Sc	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-1.76
Tb	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.16
Tm	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.10
Yb	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-1.58
Y	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-2.06

Remarks

For preparation of the standard solutions the oxides of the rare-earth elements are first dissolved in aqua regia, dried in a vacuum desiccator, and then dissolved in DMF.

With the exception of Ce, Sc, and Yb, the half-wave potentials of all metals become more positive with increasing concentrations.

³ The reference electrode is filled with supporting electrolyte.

Part 4: Actinides and technetium

Electrolytes

c(LiClO ₄) = 0.1 mol/L, pH 1.9 ... 3.1	c(LiClO ₄) = 0.1 mol/L in water. The pH value is adjusted to 1.9 ... 3.1.
c(KCl) = 0.5 mol/L, pH 2	c(KCl) = 0.5 mol/L in water. The pH value is adjusted to 2.

Half-wave potentials

Metal Ion	Supporting electrolyte	Half-wave potential	Literature
Ac	c(LiClO ₄) = 0.1 mol/L, pH 1.9 ... 3.1	-1.61	[1]
Am			[2]
Cf			[2]
Cm			[2]
Fm			[3]
Np			[4]
Pu	Al ₂ (SO ₄) ₃ / citric acid, pH 4.5	-0.43	[5]
Tc	c(KCl) = 0.5 mol/L, pH 2	-0.15	[6]

Remarks

Since all these elements are highly radioactive, it was not possible to determine them in our laboratory. The following literature references are useful.

Literature

- [1] F. David
Polarographie von Actinium
C.R. Acad. Sci. (Paris) Ser. C 271 (1970), 440-442
Ref.: Z.Anal.Chem. 256 (1971), 307
- [2] F. David
Polarography of radioisotopes 115Eu, 169Yb, 153Gd, 241Am, 244Cm, 252Cf
Compt. Rend. 270 C (1970), 2112-2115
Ref.: Electroanal. Abstr. 9 (1971), 140
- [3] F. David; M. Nussonois
Reduction of Fm(II) in acetic solution by electrolysis at controlled electric tension
Radiochem. Radioanal. Lett. 11 (1972), 1-5
Ref.: Electroanal. Abstr. 11 (1973), 328
- [4] V.F. Peretruhin; D.P. Alekseeva
Polarographic behaviour of Np(VII) in alkaline aqueous solutions
Izv. Akad. Nauk. SSSR, Ser. Khim. 23 (1974), 1438
Ref.: Electroanal. Abstr. 13 (1975), 45
- [5] G.W.C. Milner; A.J. Wood
The determination of plutonium by A.C. polarography
J. Electroanal. Chem. 7 (1964), 190-205
Ref.: Z. Anal. Chem. 219 (1965), 70
- [6] G.B.S. Salaria; C.L. Rulfs; P.J. Elving
Polarographic and coulometric determination of technetium
Anal. Chem. 35 (1963), 979-982
Ref.: Z. Anal. Chem. 219 (1965), 440

Part 5: Miscellaneous

Electrolytes

c(TBAI) = 0.1 mol/L in 80% ethanol	c(tetrabutylammonium iodide) = 0.1 mol/L in water:ethanol 80:20
c(HCl) = 0.1 mol/L	c(HCl) = 0.1 mol/L in water.
c(KSCN) = 0.1 mol/L, pH 3	c(KSCN) = 0.1 mol/L in water. The pH value is adjusted to 3
c(NaOH) = 2 mol/L	c(NaOH) = 2 mol/L in water.
HCl conc.	w(HCl) = 32%
Ca(OH) ₂ sat.	Saturated solution of Ca(OH) ₂ in water.
c(KNO ₃) = 0.3 mol/L	c(KNO ₃) = 0.3 mol/L in water.
c(H ₂ SO ₄) = 0.1 mol/L	c(H ₂ SO ₄) = 0.1 mol/L in water.
c(Na ₂ -NTA) = 0.1 mol/L	c(nitritotriacetate disodium salt) = 0.1 mol/L in water.
c(NaOH) = 0.1 mol/L / c(KNO ₃) = 0.1 mol/L	c(NaOH) = 0.1 mol/L + c(KNO ₃) = 0.1 mol/L in water.
c(EDTA) = 0.1 mol/L	c(ethylenediaminetetraacetic acid) = 0.1 mol/L in water.
c(TEA-CIO ₄) = 0.1 mol/L in DMF	c(tetraethylammonium perchlorate) = 0.1 mol/L in dimethylformamide.

Metal Ion	Supporting electrolyte	Half-wave potential
Al(III)	c(TBAI) = 0.1 mol/L in 80% ethanol	-1.67
Au(III)	c(HCl) = 0.1 mol/L	+0.90 ⁴
Ga(III)	c(KSCN) = 0.1 mol/L, pH 3	-0.84
Hf		NR
Ir(IV)	c(NaOH) = 2 mol/L	+0.38 ⁵
Nb(V)	HCl conc.	-0.60
Os(VIII)	Ca(OH) ₂ sat.	-0.42
Pt(IV)	c(KNO ₃) = 0.3 mol/L	-0.93 ⁶
Re(VII)	c(H ₂ SO ₄) = 0.1 mol/L	-0.66
Rh(III)	c(Na ₂ -NTA) = 0.1 mol/L	-0.46
Ru(I)	c(NaOH) = 0.1 mol/L / c(KNO ₃) = 0.1 mol/L	-0.81 w ⁷
Ta		NR
Ti(IV)	c(EDTA) = 0.1 mol/L	-0.34
Th(III)	c(TEA-CIO ₄) = 0.1 mol/L in DMF	-1.46 -1.71

⁴ With glassy carbon working electrode.

⁵ With platinum working electrode.

⁶ Negative peaks.

⁷ With respect to the Hg/Hg₂Cl₂/KNO₃ sat. reference electrode.