

Application Bulletin



Of interest for:
Pharmaceuticals,
food analysis

No. 213/1 e


Polarographic determination of nicotinamide

Summary

- ▶ This work is one of a series of Application Bulletins published at varying intervals for the polarographic determination of vitamins. The determination of nicotinamide, a vitamin of the B series, is described. Instructions for the determination in solutions (e.g. fruit juice), vitamin capsules and multivitamin tablets are given. The linearity range of the determination is also specified. The limit of detection is ca. 1 µg nicotinamide / 20 mL cell volume.

Theory

- ▶ Following Jacobsen and Thorgesen [3], nicotinamide is reduced at the DME according to the following scheme:



Apparatus

- ▶ 2.646.003X VA Processor with 2.647.0020 VA Stand or 2.506.0010 Polarecord or 2.626.0010 Polarecord with 2.663.002X VA-Stand

Reagents

(Except for the NaOH, these should be freshly prepared daily.)

- ▶ **Primary solution:** 1% v/v tetramethylammonium hydroxide (TMAOH). A 25% TMAOH solution (8 mL) (or 20 mL of a 10% solution) is pipetted into a 200 mL volumetric flask and this is filled to the mark with dist. water.
- ▶ **Extraction solution:** $c(\text{NaOH}) = 0.1 \text{ mol/L}$
4g/L NaOH in H₂O dist.
- ▶ **Standard, stock solution:** Ca. 2 g nicotinamide are dried overnight in a desiccator. A sample (250.0 mg) is weighed into a 250 mL volumetric flask, dissolved in dist. water and the flask filled to the mark.
- ▶ **Standard, working soln:** This is prepared as needed from the stock solution by dilution with dist. water.

Sample preparation

- ▶ **Solutions**
Solutions (injection solutions, fruit juices, etc.) are adjusted to a pH of ca. 10...12 with NaOH and filtered if need be.
- ▶ **Vitamin capsules**
The contents of 10 ... 20 capsules are weighed and mixed. The average weight is calculated and this (average capsule content) is weighed into a beaker. Extraction solution (20 mL) is added and extraction performed for 20 min on a water bath at 40°C with stirring. After cooling, the mixture is rinsed into a 250 mL volumetric flask with dist. water, the flask is filled to the mark and the contents thoroughly mixed. The flask contents are then allowed to settle for at least 20 min.

Polarographic determination of nicotinamide

Sample preparation (continued)	<ul style="list-style-type: none"> ▶ Vitamin tablets (multivitamin tablets) 10 tablets are weighed to obtain the average weight and then pulverized (grinder, mortar, etc.). The average weight of a tablet is weighed into a beaker and the procedure described for vitamin capsules then followed.
Method	<ul style="list-style-type: none"> ▶ Sample solution (1.5 ... 10 mL), which can contain 50 ... 500 µg nicotinamide, is pipetted into a polarographic vessel, 10 mL primary solution added and the solution made up to 20 mL with dist. water if need be. ▶ After deaeration with nitrogen, the DP polarogram is recorded (amplitude ca. -40 mV to -50 mV) at the DME between -1.4 V and -1.9 V. The peak potential of nicotinamide is at ca. -1.70 V. ▶ The content is determined by duplicate standard addition. It must be ensured that the concentration of nicotinamide in the cell, including that due to the standard additions, does not exceed the linearity range. ▶ Figs. 1 to 5 show the parameter settings at the 646 VA Processor, Fig. 6 and 7 an example of a nicotinamide determination.
Linearity	<ul style="list-style-type: none"> ▶ Fig. 8 and 9 show graphs nA/µg nicotinamide. They are intended to illustrate the approximate slope of the curves and should not be regarded as calibration curves.
Remarks	<ul style="list-style-type: none"> ▶ Optimum values are obtained when 50 ... 500 µg nicotinamide are present in the polarographic vessel (can be adjusted using the sample size). ▶ It is essential to use the DME as working electrode. With the SMDE lower sensitivities are obtained. The SMDE is also more susceptible to disturbances in the alkaline electrolytes (pH ca. 12.8) used.
Literature	<ul style="list-style-type: none"> ▶ Dewjatnin, W.A. / Kuznetzowa, L.A. <i>Polarographic determination of vitamins B-1, B-2 and nicotinamide in mono- and polyvitamin preparations.</i> Med. Prom. USSR 58, (1964) 58-60 (in Russian) Ref: Electroanal. Abstr. 5, (1967) 43 ▶ Göbbeler, K.H. / Breinlich, J. <i>Quantitative, wechselstrompolarographische Simultanbestimmung von Vitaminen der B-Gruppe.</i> Pharm. Ztg. 48, (1972) 1859-1862 ▶ Jacobsen, E. / Thorgersen, K.B. <i>Electroreduction and pulse-polarographic determination of nicotinamide in multivitamin tablets.</i> Anal. Chim. Acta 71, (1974) 175-184 ▶ Moore, J.M. <i>Polarographische Schnellanalyse von Nikotinamid in pharmazeutischen Präparaten.</i> J. Pharm. Sci. 58, (1969) 1117-1120 (in English) Ref: Fresenius, Z. Anal. Chem. 254, (1971) 159 ▶ Söderhjelm, P. / Lindquist, J. <i>Electrochemical assay of thiamine, riboflavine, pyridoxine, nicotinamide and ascorbic acid in pharmaceutical preparations.</i> Acta Pharm. Suec. 13, (1976) 201-212 ▶ Taira, A.Y. <i>Polarographic determination of niacinamide in multivitamin preparations.</i> J. Assoc. Off. Anal. Chem. 57/4, (1974) 910-913 Ref: Electroanal. Abstr. 13, (1975) 388

Polarographic determination of nicotinamide

Fig. 1 Example program page 2 646 VA Processor

1 Determination of niacinamide in unknowns	METHOD 4 PAGE 2
2 MPL 1 EL.TYPE MME	GEN.SPECIFICATIONS
PARAMETERS	
3 iR.MODE	N
4 SPEED	5
5 D.SIZE	5
6 N.DROPS	5
RECOGNITION	
7 SPIKE THRESH	5
8 H.THRESH	4
9 U.TOL	7
10 W.TOL	8
11 ASYM.TOL	9

Fig. 2 Example program page 3 646 VA Processor

Determination of niacinamide in unknowns	METHOD 4 PAGE 3
MPL 1 EL.TYPE MME	OPERATION SEQUENCE
OPERATIONS/PARAMETERS	
1 PURGE ;STIR ;	5 s
2 [ADDL ;OPURGE;OSTIR ;	5 s
3 (REP ;	
4 DME ;MEAS ;	5 s
4a M.MODE DPN	-40 mV
4b T.STEP	600 ms
4c U.SET	-1.400 V
5 SWP 0 ;	50 s
5a U.END	-1.900 V
5b U.STEP	6 mV
SW.RATE	10.0 mV/ s
6 OMEAS ;	
7 REP) 1;	
8 OMEAS ;PURGE ;STIR ;	
9 BEEP ;ADD1]2;	30 s
10 OMEAS ;OPURGE;OSTIR ;	
11 BEEP ;END ;	

Fig. 3 Example program page 4 646 VA Processor

Determination of niacinamide in unknowns	METHOD 4 PAGE 4				
MPL 1 EL.TYPE MME	ALLOCATIONS				
a SOLUTE	b U.VERIF	c DOS	d V.SOLN	e m.CONC	f m.BLANK
Subst	Ux	SoIn	c. v	rho.x	bx
1 Niamde	-1.700 V	1	c 100 uL	1.000 g/ L	0.000 g
2					
3					
4					
5					
6					
7					
8					
9 SUPP.ELEC	1%TMAH+DDI pH 13.0				
10 V.MEAS	20.000 mL				
11 ALIQUOT	1.000				
12 DATE	91-04-02				
13 TIME	10:10				

Polarographic determination of nicotinamide

Fig. 4 Example program page 5 646 VA Processor

```

Determination of niacinamide in unknowns
MPL 1          EL.TYPE MME
              a          b          c          d          e
SEGMENT        Y.AXIS/L  Y.AXIS/R  X.AXIS/DIV
1 SWP 0        0.00 A    -1.00 uA   100 mV
2
3
4
5
6
7
8

9 RECORD      SXXX      FR
10
11 SEND
12
  
```

Fig. 5 Example program page 6 646 VA Processor

```

Determination of niacinamide in unknowns
MPL 1          EL.TYPE MME
              a          b          c          d
ANALYTE        EVAL      R.QUANT  R.UNIT      SIGNIF.DIG
1 Niamde       N         rho(Niamd) mg/g        4
2
3
4
5
6
7
8

          (EV.QUANT  + ADDEND)  * FACTOR  / DIVISOR
11 Niamde   A         0.00000    532.420   1.00000
12
13
14
15
16
17
18
  
```

Polarographic determination of nicotinamide

Fig. 6 Curve examples: Determination of nicotinamide in tablets of declared content 100 mg/tablet




Fig. 7 Result block of the determination from Fig. 6

```

METROHM 646 VA-PROCESSOR (5.646.8041)
Determination of niacinamide in unknowns      METHOD 4
MPL 1      EL.TYPE    MME


SUPP.ELEC    1%TMAH+DDI pH 13.0
V.MEAS      20.000 mL
ALIQOT      1.000

REMARK      Detrn. of niacinamide in niacinamide tablets
            Ag/AgCl (3M KCl) reference electrode
NAME      Prof.J.G.Dick
RUN#      2

ANALYTE    L R S      U.SUBST    EV.VALUE    DELTA      m.ANALYTE
Niamde    A0 0 0    -1.707 V    373.6 nA
          A0 1 0    -1.708 V    376.1 nA
          A1 0 0    -1.708 V    587.0 nA
          A1 1 0    -1.709 V    587.9 nA    212.6 nA
          A2 0 0    -1.709 V    792.1 nA
          A2 1 0    -1.709 V    794.7 nA    205.9 nA
          m.STD    100.0 ug    SLOPE    477.7 ug/uA    179.6 ug

rho(Niamd) =    95.64      mg/g

SMPL.V,m      1.00000 mg      IDENT Swiss Herbal vit.pp
DATE 91-04-02    TIME 09:31
    
```

Polarographic determination of nicotinamide**Fig. 8** Linearity curve 25 ... 225 μg nicotinamide / 20 mL**Fig. 9** Linearity curve 100 ... 900 μg nicotinamide / 20 mL